

SPIDER SCREW[®]



Temporary orthodontic anchorage system



WYŁĄCZNY DYSTRYBUTOR W POLSCE

Ortotech
Orthodontic Materials
and Technologies

ZAMÓWIENIA:

tel: 22 10 07 207

mob: 889 779 448

e-mail: info@ortotech.pl

web: www.ortotech.pl

ul. Juliana Smulikowskiego 10/13
00-389 Warszawa

SPIDER SCREW[®]



Temporary orthodontic anchorage system

A close-up photograph of a young woman with blonde hair, smiling broadly, showing her teeth. The image is partially obscured by a blue geometric shape that frames the top and right sides of the page.

Nothing has been left to chance.


The Spider Screw has obtained an **international patent** since its inception, due to its innovative characteristics: the simultaneous presence of external and internal rectangular slots and round internal slots.

The **SPIDER SCREW'S** geometry is a result of careful design in each of its three components: the **ORTHODONTIC HEAD**, the **TRANSMUCOSAL PORTION** and the **INFRABONY PORTION**.

The **SPIDER SCREW** has been developed to offer a number of versatile anchorage options capable of immediate loading. Immediate loading is possible because the **SPIDER SCREW** is a non osteointegratable implant and consequently force can be applied immediately after placement. The applied force can range from 50 to 300 grams depending on screw choice, bone quality, and the desired orthodontic movement.

The **SPIDER SCREW'S** design is extremely versatile, due to its small dimensions and unique design. It is easily placed in either the maxilla or mandible, even where access is limited and bone quality is less than ideal. Placement is simplified by the self-drilling feature found in the K1 and K2 systems.

The **SPIDER SCREW** is an anchorage device that can be used during every phase of orthodontic treatment and is suitable for symmetric or asymmetric anchorage. The screws assists in the success of orthodontic treatment, both in adults and adolescents, by reducing treatment time without patient co-operation.



The **SPIDER SCREWS** are provided in different diameters and lengths so to be placed in different sites (edentulous ridges, palate, tuberosity, interdicular areas, infrazygomatic crest, retromolar region, mandibular ramus and synfisis).

The **SPIDER SCREW** is made in Grade 5 Titanium alloy, supplied clean and sterile inside a sealed polypropylene bottle. The screws have been cleaned, decontaminated and sterilized via gamma radiation. Included in the package are three removable stickers containing important information (device name, reference code, lot number, etc.) which can be applied to the patient's record card for traceability.



1. **ORTHODONTIC HEAD** (fig. 1,2)

The orthodontic head was designed to facilitate appliance (wires, springs, etc.) placement. There is a **bracket like head featuring two intersecting .022" slots**. The under tie wing area can also function as another .022" x .025" slot and features **two intersecting slots of .027"** in diameter with chamfered inlets to simplify insertion of wires or ligatures. The head is available different dimensions: the largers (Spider Screw K2 Regular Plus 2 mm, Spider Screw K2 1,9 mm and Spider Screw Konic 2/2,3 mm) are for palatal use; the smaller (Spider Screw K1 1,5 mm, Spider Screw C1 1,5 mm, and Spider Pin 1,3 mm) are for patient comfort.



2. **TRANSMUCOSAL PORTION** (fig. 3,4)

The length of the transmucosal portion is variable and allows for optimal adaptation to different intraoral mucosa thicknesses during bio-maintenance. **Short** for areas of thin attached gingiva. **Long** for areas with thick or freely moveable tissues. The transmucosal portion is polished with a special treatment to help avoid soft tissue irritation and make cleaning easier to accomplish.



3. **INTRABONY PORTION** (fig. 5)

The Spider Screw's thread shape has an asymmetrical profile making it easy to place while ensuring maximum stability and avoiding bone stress. Spider Screw K1 1,5 mm, Spider Screw K2 1,9 mm and Spider Screw Self-Ligating, tapered thread, are self-drilling and self-tapping which makes pre-drilling before insertion unnecessary - dependant upon bone structure. This makes the **Spider Screw K1, K2** and **Self-Ligating** easy to place while reducing the risk of root damage. Spider Screw C1 1,5 mm and Spider PIN 1,3 mm cylindrical thread, require pre-drilling and are used in areas that have poor bone quality or greater retention requirements.



Insertion Sites:

MAXILLA

- Infrazygomatic crest
- Edentulous ridges
- Palate
- Tuberosity
- Interadicular areas

MANDIBULA

- Edentulous ridges
- Retromolar region
- Mandibular ramus
- Interadicular areas
- Synfisis

INDICATIONS

SPIDER SCREW ANCHORAGE SYSTEM (SSAS) allows sagittal and vertical movement of all teeth (intrusion, extrusion, distalization and mesialization) and can be used for treating the following:

- > Class I, Class II, Class III malocclusion treatment;
- > Anchorage Recovery;
- > Anchorage Reinforcement;
- > Asymmetrical case management;
- > Uprighting of upper and lower molars;
- > Correction of over erupted teeth (molar, premolar, incisors);
- > Deep bite and open bite conditions;
- > Pre-Prosthetic Orthodontic treatment;
- > Borderline cases;
- > Orthodontic treatment without patient cooperation (MBGM system).

GENERAL INFORMATION

The placement of **SPIDER SCREW** is a procedure requiring specific knowledge of anatomy and technique. It is absolutely necessary that it is carried out by specifically trained doctors. It is important to know that improper patient selection and/or incorrect technique can cause placement failure and/or loss of supporting bone.

An effective and complete screening of the patient must be performed and each case carefully evaluated. A very thorough examination is needed, as well as anatomical reference for the evaluation of bone quantity and quality using radiographic research (Long Cone Endoral Radiograph, Orthopantograph, Teleradiography, and Computerized Tomography).

Carefully read the instructions for use inside the package before the Spider Screw placement. The **SPIDER SCREW** is for single use only and should not be reused. Use only the instruments mentioned in this catalog, making sure that all the instruments are sterilized and efficient.

It is suggested to disinfect the insertion area and give local anesthesia as needed. It is very important that the clinician attends a training course for a complete overview of all the possible applications, as this catalog shows only a few.

TAPERED THREAD AND SELF-LIGATING SCREW

If a **SPIDER SCREW** is to be inserted in an edentulous area where there is bone availability, references from a panoramic radiograph can be sufficient.


1. In areas close to delicate anatomical structures, such as interadicular spaces, a long cone radiograph is recommended.
2. A surgical splint can be made with orthodontic wire, fixing it to the teeth with acrylic or thermoplastic resin. The orthodontic wire is inserted in the acrylic resin and is appropriately bent so that its tip corresponds to the point of insertion of the **SPIDER SCREW**.
3. Use a periapical radiograph (by using the long-cone parallel technique) to verify the correct placement of the orthodontic wire.
4. The insertion site can be marked with a pressure point or methylene blue dot on the soft tissue. In mobile mucosa it is recommended to leave the surgical guide in place during the drilling phase and/or the screw insertion.



1.




2.



3.



4.



5. After site disinfection (chlorhexidine) insert the **SPIDER K1** or **K2** using the **pick-up driver shaft DSP-5052S** (insert the Spider K1 or K2 **SELF-LIGATING** using the **pick-up driver shaft DSP-5652S**). It is also possible to use the **contra-angle pick-up driver at low speed 25/30 rpm** (**DPQ-2820** for the Spider Screw K1 or K2 Tapered Thread, **DPQ-3420** for the Spider Screw Self-Ligating).

In order to avoid excessive torque stress during insertion, which could cause bone compression and consequent recession or cause the screw to break, it is recommended to use a technique of alternating between screwing and unscrewing to gradually ease the screw into position. Final placement is achieved by using the handle driver **DSQ-2824** for the Spider Screw K1 or K2 Tapered Thread or **DSQ-3424** for the Spider Screw Self-Ligating, to complete the insertion as this provides the most controlled tactile method.

6. In the case of very compact bone use a spiral drill (**FSC-1108** for K1 Tapered Thread and Self-ligating or **FSC-1309** for K2 Tapered Thread and Self-ligating) to make a pilot hole which makes screw insertion easy to perform.

SPIDER SCREW C1 AND PIN PLACEMENT

Follow points 1 to 4 as above.

- 5a. After site disinfection (chlorhexidine) the spiral drill is used to perforate the soft tissue and cortical bone (no incision needed). Cold irrigation is used during the drilling procedure (5°C/41°F). Use the \varnothing 0,9 mm drill for the **SPIDER SCREW PIN** and the \varnothing 1,2 mm drill for the **C1**.

- 6a. You can choose between two options: manual or mechanical insertion. For manual insertion use **DSP-2352S** for the Spider Screw PIN \varnothing 1,3 mm and **DSP-5052S** for the Spider Screw C1 \varnothing 1,5 mm. For mechanical insertion use the contra angle pick-up driver **DPQ-2352S** for the Spider Screw PIN \varnothing 1,3 mm and **DPQ-2820** for the Spider Screw C1 \varnothing 1,5 mm mounted on a low speed contra-angle handpiece (25/30 rpm). Final placement is achieved by using the handle driver to complete the insertion as this provides the most controlled tactile method.

POST APPLICATION PATIENT INSTRUCTION

Application of chlorhexidine rinse 2 – 3 times per day for the first 7 days.

Perform normal hygiene procedures. The patient should brush the screw normally as a tooth.

SPIDER SCREW REMOVAL

To remove the **SPIDER SCREW**, it is simply unscrewed with the appropriate screwdriver. For anterior and lateral areas is advisable to use the **handle driver**. While for posterior areas is advisable to unscrew with **contra angle pick-up driver**. This can be accomplished with or without anesthesia.

If the Spider Screw does not unscrew easily it is recommended to use a technique of alternating between unscrewing and screwing. Healing takes place in a few days.



SPIDER SELF-LIGATING K1

Ø1,5 mm


Tapered Thread (Self-tapping / Self-drilling)

Available in length 6,5 – 8 – 10 mm.

The **SPIDER SCREW SELF-LIGATING TAD - K1** is self-drilling and self-tapping. Due to the design of the tapered thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the Ø 1,1 mm drill provided to penetrate the cortical plate.

During insertion of the Spider K1 do not exceed 20 N/cm.

SXL-1506	Ø 1.5 x 6.5 mm
SXL-1508	Ø 1.5 x 8.0 mm
SXL-1510	Ø 1.5 x 10.0 mm



SPIDER SELF-LIGATING K2

Ø1,9 mm


Tapered Thread (Self-tapping / Self-drilling)

Available in length 6 - 7 - 9 - 11 mm.

The **SPIDER SCREW SELF-LIGATING TAD - K2** is self-drilling and self-tapping. Due to the design of the tapered thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the Ø 1,3 mm drill provided to penetrate the cortical plate.

During insertion of the Spider K2 do not exceed 30 N/cm.

SXL-1906	Ø 1,9 x 6mm
SXL-1907	Ø 1.9 x 7 mm
SXL-1909	Ø 1.9 x 9 mm
SXL-1911	Ø 1.9 x 11 mm



CSS-6008
 SELF LIGATING
 SPIDER SCREW ORGANIZER



FSC-1108 DRILL Ø 1,1 mm - K1



FSC-1309 DRILL Ø 1,3 mm - K2



DPQ-3420 CONTRA ANGLE PICK-UP DRIVER - SHORT



DPQ-3425 CONTRA ANGLE PICK-UP DRIVER - LONG



DLM-3134 DRIVER
 for closing the slot



DXL-2820 SELF LIGATING TAD KEY



DSQ-3424 HANDLE SQUARE DRIVER



DSP-5652S PICK-UP DRIVER SHAFT



DSX-1690S SCREWDRIVER BODY

SPIDER LINK

A system of plates for skeletal anchorage


- > Allows application of mini-screw in more suitable anatomical areas.
- > Easy to place and remove.
- > Easy to mould.

SPIDER LINK is a skeletal anchorage system comprising of a SELF-LIGATING SCREW and titanium preformed orthodontic anchorage attachment.


Current skeletal anchorage systems don't always allow the insertion of screws in the desired position, due to teeth interference and/or anatomy.

This new system, SPIDER LINK (SELF-LIGATING SPIDER SCREW + ORTHODONTIC ANCHORAGE ATTACHMENT) lets you overcome this limitation.


The screw can be placed in more suitable anatomical areas and then, thanks to the orthodontic anchorage attachment, it is possible to apply forces to desired teeth.




PSLV-0003
V-FORM




PSLV-0004
V-FORM




PSLH-0005
H-FORM



PSLH-0001
H-FORM




PSLH-0003
H-FORM




PSLH-0004
H-FORM


CLINICAL CASES




Orthodontic Anchorage Attachment for molar Mesialization and Distalisation.




Orthodontic Anchorage Attachment for molar Distalisation.




Orthodontic Anchorage Attachment for open bite condition.
Second molars and molars intrusion.



Orthodontic Anchorage Attachment for open bite condition.
Asymmetrical intrusion.



Class II Division 1 crowding treatment
Orthodontic Anchorage Attachment for simultaneous bilateral molar distalisation.



Class II Division 1 on the left side.
Orthodontic Anchorage Attachment for simultaneous molar distalisation.



SPIDER K1

Ø1,5 mm

Tapered Thread (Self-drilling/Self-tapping)

Available in length 6,5 – 8 – 10 mm.

The **SPIDER SCREW K1** is self-drilling, due to the design of the tapered thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the Ø 1,1 mm drill provided to penetrate the cortical plate. Available in the following versions:

LONG NECK - elongated neck height (2 mm) for areas of thick tissue (Posterior and lateral areas).

SHORT NECK - standard neck height (1 mm) for areas of thin tissue (Anterior and lateral areas).

	SCL-1506	Long Neck Ø 1.5 x 6.5 mm
	SCL-1508	Long Neck Ø 1.5 x 8.0 mm
	SCL-1510	Long Neck Ø 1.5 x 10.0 mm
	SCR-1506	Short Neck Ø 1.5 x 6.5 mm
	SCR-1508	Short Neck Ø 1.5 x 8.0 mm
	SCR-1510	Short Neck Ø 1.5 x 10.0 mm



SPIDER K2

Ø1,9 mm

Tapered Thread (Self-drilling/Self-tapping)

Available in length 6 - 7 - 9 - 11 mm.

The **SPIDER SCREW K2** is self-drilling, due to the design of the tapered thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the Ø 1,3 mm drill provided to penetrate the cortical plate. The 5 mm and 6 mm lengths screws can be placed in the palate or in areas where the bone thickness is reduced. Available in the following versions:

LONG NECK - elongated neck height (2 mm) for areas of thick tissue (Posterior and lateral areas).

SHORT NECK - standard neck height (1 mm) for areas of thin tissue (Anterior and lateral areas).

	SCL-1906	Long Neck Ø 1.9 x 6.0 mm
	SCL-1907	Long Neck Ø 1.9 x 7.0 mm
	SCL-1909	Long Neck Ø 1.9 x 9.0 mm
	SCL-1911	Long Neck Ø 1.9 x 11.0 mm



SCR-1906	Short Neck Ø 1.9 x 6.0 mm
SCR-1907	Short Neck Ø 1.9 x 7.0 mm
SCR-1909	Short Neck Ø 1.9 x 9.0 mm
SCR-1911	Short Neck Ø 1.9 x 11.0 mm



SPIDER C1

Ø1,5 mm

Cylindrical Thread (Pre-drilling/Self-tapping)

Available in length 6,5 – 8 – 10 mm.

The **SPIDER SCREW C1** is self-tapping and requires pre-drilling with a drill of diameter 1,2 mm. Available in the following versions:

LONG NECK - elongated neck height (2 mm) for areas of thick tissue (Posterior and lateral areas).

SHORT NECK - standard neck height (1 mm) for areas of thin tissue (Anterior and lateral areas).



SLP-1506	Long Neck Ø 1.5 x 6.5 mm
SLP-1508	Long Neck Ø 1.5 x 8.0 mm
SLP-1510	Long Neck Ø 1.5 x 10.0 mm



SSM-1506	Short Neck Ø 1.5 x 6.5 mm
SSM-1508	Short Neck Ø 1.5 x 8.0 mm
SSM-1510	Short Neck Ø 1.5 x 10.0 mm



CSS-4009
SPIDER SCREW ORGANIZER



FSC-1210

DRILL Ø 1.2 mm - C1



FSC-1108

DRILL Ø 1.1 mm - K1



FSC-1309

DRILL Ø1.3 mm - K2



DPQ-2820

CONTRA ANGLE PICK-UP DRIVER – SHORT



DPQ-2825

CONTRA ANGLE PICK-UP DRIVER – LONG



DPX-2830

CONTRA ANGLE CROSS DRIVER



DSQ-2824

HANDLE DRIVER SHORT



DPH-2824

PICK-UP HANDLE DRIVER



DSX-2852S

CROSS DRIVER SHAFT
for DSX-1690S



DSP-5052S

PICK-UP DRIVER SHAFT
for DSX-1690S



DSX-1690S

SCREWDRIVER BODY


CLINICAL CASES



Intrusion Posterior Areas



Direct Anchorage – Molar Uprighting and Intrusion



Indirect Anchorage (by kind permission of Dr. B.G. Maino)



Molar Uprighting and Intrusion (by kind permission of Dr. N.Derton)



SPIDER PIN

Ø 1,3 mm

Cylindrical thread (Self-Tapping)

Available in length 8 – 10 mm.

SPIDER PIN is a simple design Self-Tapping screw and requires pre-drilling. It is ideal for areas where a reduced size head is required (i.e. narrow interproximal spaces). The Spider PIN head is rounded to increase patient comfort and allow for better cleaning.

- Ideal for narrow interproximal spaces.
- Simple head, perfect for NiTi closed coil spring attachments or elastic chains.
- Requires no patient cooperation and reduces treatment time.
- Smooth, rounded design for patient comfort.

SCL-1308	Ø 1.3 x 8.0 mm
SCL-1310	Ø 1.3 x 10.0 mm

CLINICAL CASE



Spider Pin Application

CSS-3006
SPIDER PIN ORGANIZER





FSC-0910

DRILL - PIN



DPQ-2322

CONTRA ANGLE PICK-UP DRIVER



DSQ-2324

HANDLE DRIVER – SHORT



DSP-2352S

PICK-UP DRIVER SHAFT
for DSX-1690S or DST-1600



DSX-1690S

SCREWDRIVER BODY

OPTIONAL DEVICES WORKING WITH ALL SYSTEMS



CAM-2000

MANUAL CONTRANGLE



DST-1600

TORQUE DRIVER

- Torque driver to control insertion force.
- The torque driver has a range of torque of 15 N/cm: from 5 up to 20 N/cm.
- Interchangeable shafts for all kind of Spider Screw.
- Can be sterilized.

SUMODIS

Simultaneous upper Molars Distalizing System

SUMODIS it's a system for simultaneous distalization of upper molars in presence of second molars in II class non-extractive treatments without co-operation.

SUMODIS is a combination of sliding mechanics using Spider Screws as a sole anchorage resource.

SUMODIS (Simultaneous Upper Molar Distalizing System) main characteristic is two distalizing components: one against the first molar and the other one against the second, working simultaneously.

One Kit Sumodis contains:

Stainless steel tempered wire .016x.022

Neosentalloy wire .018x.025 provided with stop

Protection elastomer tube

Neosentalloy open coil spring 200 g.

Crimpable stop

Double tube


Preshaped palatal bar Ø 1 mm




CLINICAL CASE




Sumodis in place




First distalizing component. Insert in a sectional .016" x .022" stainless steel wire and a 200 g. NeoSentalloy open coil between the first premolar and the molar.




Before fixing .016" x .022" stainless steel wire, insert the double tube into it using the lower tube and position it near by the premolar bracket. The double tube, which can slide, is blocked by the compressed coil on one side and the premolar bracket on the other side.




The second distalizing component is a shape memory .018" x .025" stainless steel NeoSentalloy wire (with crimped stop). Put the second stop at a distance of 6 mm longer than the distance between the mesial side of the buccal tube on the second molar and the distal side of the double tube on the sectional wire.



When the NeoSentalloy wire is inserted into the tube on the second molar and in the double tube free hole, it will raise in the buccal fold and activate 6 mm.



In this phase you can see both distalizing components simultaneously activated.



End of first phase MBGM system. Molars Distalization done.



SPIDER K2 REGULAR PLUS

Ø 2 mm

Cylindrical Thread (Self-drilling)

Available in length 6 – 7 – 9 – 11 – 13 – 15 mm.


The inner threaded portion of the **SPIDER SCREW K2 REGULAR PLUS** opens up vast possibilities for applications.

Additional devices can be fitted to the screw, such as acrylic resin abutments, stainless steel abutments, pre-shaped abutments or lab made units.


The palatal insertion of the **SPIDER SCREW K2 REGULAR PLUS** and the possibility of preparing specific orthodontic devices attachable to it allows an innovative way of treating class II and class III conditions, reducing treatment time and requiring no patient co-operation.

SSP-2006	Regular Plus Ø 2 x 6.0 mm
SSP-2007	Regular Plus Ø 2 x 7.0 mm
SSP-2009	Regular Plus Ø 2 x 9.0 mm
SSP-2011	Regular Plus Ø 2 x 11.0 mm
SSP-2013	Regular Plus Ø 2 x 13.0 mm
SSP-2015	Regular Plus Ø 2 x 15.0 mm


CLINICAL CASES



4.



5.



By kind permission of Dr. B.G. Maino



FSC-1319

LONG DRILL Ø 1,3 MM - RP



PMA-5008

ACRYLIC RESIN ABUTMENT



PMA-4050R

STAINLESS STEEL ABUTMENT



PMA-5006R

UTILITY ABUTMENT



CTS-5000

COPING TRANSER



ANR-4012

LAB ANALOG



DSE-2410

DRIVER



DSX-1690S

SCREWDRIVER BODY

SPIDER SCREW


Regular Plus KONIC

ORTHODONTIC HEAD



The Cylindrical head is designed to overcome possible disparallelisms between the screw inclination and the customized appliance to be used.

TRANSMUCOSAL PORTION



The transmucosal portion is polished with a special treatment for an excellent bonding to the oral mucosa and to help avoid soft tissue irritation.

THREADING



The asymmetric thread profile is specifically studied to provide an optimal screw stability and an easier screw insertion.

TIP



Conic Flute.
The Spider Screw's tip shape makes it easy to place it and to maintain the inclination angle required.

SPIDER SCREW

Regular Plus KONIC



FSC-1309G DRILL FOR SSP Ø 2,0 mm

FSC-1609L DRILL FOR SSP Ø 2,3 mm



DPQ-3825 PICK-UP DRIVER FOR GUIDED SURGERY



PRL-3513 DRILL EXTENSION



DPM-3075 CONTRA ANGLE PICK-UP DRIVER



BGU-6050 SURGICAL GUIDETUBE



TUB-1001 MESIAL TUBE WITH HOOK
TUB-1003 MESIAL TUBE WITHOUT HOOK



PSRP-1001N WIRE SUPPORT RING



PSRP-2005N PLATE WITH WIRE



VRS-1650 FIXATION SCREW FOR PLATE AND RING



CSS-7014 SPIDER SCREW Regular Plus Kit

SCREW LENGTHS Ø 2,0 mm

- SSP-2007N Regular Plus Ø 2 x 7 mm
- SSP-2008N Regular Plus Ø 2 x 8 mm
- SSP-2009N Regular Plus Ø 2 x 9 mm
- SSP-2010N Regular Plus Ø 2 x 10 mm
- SSP-2011N Regular Plus Ø 2 x 11 mm
- SSP-2012N Regular Plus Ø 2 x 12 mm
- SSP-2013N Regular Plus Ø 2 x 13 mm
- SSP-2014N Regular Plus Ø 2 x 14 mm
- SSP-2015N Regular Plus Ø 2 x 15 mm

COMPONENTS

SCREW LENGTHS Ø 2,3 mm

- SSP-2307N Regular Plus Ø 2,3 x 7 mm
- SSP-2308N Regular Plus Ø 2,3 x 8 mm
- SSP-2309N Regular Plus Ø 2,3 x 9 mm
- SSP-2310N Regular Plus Ø 2,3 x 10 mm
- SSP-2311N Regular Plus Ø 2,3 x 11 mm
- SSP-2312N Regular Plus Ø 2,3 x 12 mm
- SSP-2313N Regular Plus Ø 2,3 x 13 mm
- SSP-2314N Regular Plus Ø 2,3 x 14 mm
- SSP-2315N Regular Plus Ø 2,3 x 15 mm



STT-3016 SPRING STOP



PSRP-2001N PLATE



DSQ-2507 DRIVER



ANR-3812 REGULAR PLUS ANALOG



PMA-3850 STAINLESS STEEL ABUTMENT



CTS-6000 COPING TRANSFER



SCB-5010 SCANBODY



SRP-3850 NITI SPRING 500 GR



VRS- 1645 ABUTMENT FIXATION SCREW

By kind permission of Dr. B.G. Maino

Partner Lab: Orthomodul, CDT Emanuele Paoletto





HDC S.r.l.
 Via dei Mestieri 5/7, Thiene 36016 (VI) Italy
 Tel. +39 0445 364148
 info@hdc-italy.com

CERTIFICATIONS



HDC obtained CE marking and ISO 13485, which prove that our product has been assessed and meets EU safety, health and environmental protection requirements, as well as the requirements for a comprehensive quality management system for the design and manufacture of medical devices.

kiwa
Partner for progress

Reg. Numero: 9950-M
 Data di nascita: 2012-06-29 Data di ultima modifica: 2015-06-25
 Data di prossimo rinnovo: 2018-06-28

**Certificato del Sistema di Gestione per la qualità
 ISO 13485:2003**

Si dichiara che il sistema di gestione per la Qualità dell'Organizzazione:
H.D.C. S.r.l.
 è conforme alla norma UNI CEI EN ISO 13485:2012 per i seguenti prodotti/servizi:
 Progettazione e produzione di viti ortodontiche, impianti dentali e relativi accessori.
 Commercializzazione di prodotti per il settore ortodontico

Chief Operating Officer
 Giampiero Belcredi

Il mantenimento della certificazione è soggetto a sorveglianza annuale e subordinato al rispetto dei requisiti contrattuali Kiwa Cermet Italia.

Riferirsi al manuale qualità per i dettagli delle esclusioni ai requisiti della norma UNI CEI EN ISO 13485:2012 il presente certificato è costituito da 1 pagina.

H.D.C. S.r.l.
 Via Asiago, 44
 36030 Sarcedo VI Italia

Kiwa Cermet Italia S.p.A.
 Società con unico socio, soggetta all'articolo di direzione e coordinamento di Kiwa Italia Holding Srl
 Via Caramelo, 28
 40057 Granarolo dell'Emilia (BO)
 Tel. +39 051 438 3 111
 Fax +39 051 763 382
 E-mail: info@kiwacermet.it
 www.kiwacermet.it

CERMET

IAF

ACCREDIA 002 N° 0074 001 N° 0060
 024 N° 0160 F04 N° 0041
 PID N° 0468

**CERTIFICATO CE DEL SISTEMA DI
 GARANZIA DELLA QUALITÀ**
 EC Quality Assurance System Certificate

Si certifica che, sulla base dei risultati degli audit effettuati, il Sistema completo di garanzia di Qualità dell'Organizzazione:
 Kiwa Cermet S.p.A. on the basis of audit carried out, the full Quality Assurance System of the Company.

H.D.C. S.r.l.

Reg. No: MED 9911-A

Indirizzo / Address:
 Unità Operativa / Operational unit:
 Via Asiago, 44
 36030 Sarcedo, VI - Italia
 Sede legale / Registered Headquarter:
 Via dell'Industria, 19
 Sarcedo (VI) - VI - Italia

È conforme ai requisiti applicabili della / Is in compliance with the applicable requirements of:
**Direttiva 93/42/CEE e successive modifiche ed integrazioni, Allegato II escluso il pto 4,
 attuata in Italia con Digs. 46 del 1997/02/24 e successive modifiche ed integrazioni**
 93/42/EEC Directive as amended, Annex II without point 4, transposed in Italy by Digs. 46 of 1997/02/24 as amended
 per le seguenti tipologie di Dispositivi Medici / for the following Medical Devices:

**Componentistica protesica e accessori / Prosthetic components and accessories
 Frese e strumentario odontoiatrico / Drills and dental equipment
 Impianti dentali / Dental implants**

Identificazione / Identification: Vedere allegato tecnico al presente Certificato / See technical sheet enclosed to this certificate

Il presente Certificato è soggetto al rispetto del Regolamento CERMET ed è valido solo per le tipologie di dispositivi sopra identificate soggette a sorveglianza.
 This Certificate is subject to CERMET regulations and it is valid only for the above mentioned Medical Devices that are subject to survey.
 The technical sheet is an integral part of this Certificate.

Data di emissione / Issue date: 2015/02/20
 Data ultima modifica / Last revised date: 2015/02/20
 Data scadenza / Expiry date: 2020/02/19
 Revisione / Revision: 0
 Pagina / Page: 1 di 3

CE Organismo Notificato n. 0476
 European Notified Body n. 0476

Direttore Generale
 General Manager
 Giampiero Belcredi

kiwa

CERMET
 Certificazione e ricerca per la qualità

Kiwa Cermet Italia S.p.A. - Sede legale - Via Calabro 21 - 41057 Capriano di Carpi (MO) - Tel. +39 051 438 111 - Fax +39 051 763 382 - www.kiwacermet.it